

Reduction and Analysis graphic user interface (RAGui)

Shun Yu

Goal: development and simple graphic user interface, which can reduce the scattering data in hdf5 format, offer multiple plotting options and perform simple line profile analysis

Available Functions:


- Easy loading and viewing raw data in hdf5
- Data reduction with multiple function selection
- Saving reduction data in hdf5
- Fast 1D plot with multiple choice:
 - Single plot within one hdf5
 - Compare single line profile from different hdf5
 - Waterfall plot within one hdf5
- Fast 2D plot of the “cake” cut and profile change as a function of frame number

Under development:

- Multicomponent fitting
- Plotting as a function of meta data, e.g. motor position and time stamp

1. Define raw data path, saving path, and basic GUI operation

2. Tab menu to selection between Raw data viewing and 1D profile plotting



1. Define raw data path, saving path, and basic GUI operation



- “Raw Data”: Define the raw data path
- “Exit”: quit GUI

Default path:


Raw Data → “./test data”

Saving Data → “./processed data”


Feel free to change the path



- “Saving Path”: Define the saving data
- “Calibration”: call “pyFAI-calib2” Calibration GUI




- “Help wiki”: Open browser for User Manual
- “About”: introduction of the GUI



Mask and Poni file could be created from this GUI





1. Raw data File handling

- Refresh: user defined quick reload from Raw data folder (define in GUI Menu "File") and show in file in [Panel 2](#)
- Load: open pop up window to load any hdf5 from any accessible folder and show in file in [Panel 2](#)
- Remove: user selected single file or multiple file in [Panel 2](#) (occasionally Exception warning pop up, but does not influence the use once warning is acknowledged)
- Clear: remove all listed file from [Panel 2](#)

2. hdf5 tree view:

Visualized all the structure inside readable hdf5

3. Data view window:

Default silx viewer's data viewing window with multiple choice of showing the data

4. Reduction Panel:

Three checkable Groupbox: *Process*, *Remove Background* and *Integration*

- *Process*: Normalization function (**Selection include: transmission (I_t), flux (I_0), exposure time(Δt) and all of them**)
- *Remove Background*: select background in **hdf5 (cosaxs format)** and calculate the **average** over all the detector image as the background with user defined scaling factor of the background. Removing background is done by **Subtracting** background from selected single or multiple file in [Panel 2](#)
- selected single or multiple files *Integration*: azimuthal integration of all detector images in selected single or multiple from [Panel 2](#).

5. Combo processing



Combined Reduction to perform all Checked Groupbox functions **in Panel 4 in one click.**


(Working but not implemented at CoSAXS yet)

6. ProgressBar

It shows the progress of actions in [Panel 4](#) and [5](#) for each file (despite multiple files are selected)

Select between 1D and 2D plot





1. Reducted File handling

- Refresh: user defined quick reload from Raw data folder (define in GUI Menu "File") and show in file in **Panel 2**
- Load: open pop up window to load any hdf5 from any accessible folder and show in file in **Panel 2**
- Remove: user selected single file or multiple file in **Panel 2** (occasionally Exception warning pop up, but does not influence the use once warning is acknowledged)
- Clear: remove all listed file from **Panel 2**

2. hdf5 tree view:

Visualized all the structure inside readable hdf5

3. 1D profile plotting selection

- Selected between *Q plot* and *Azimuthal plot*

Q_{low} and Q_{high} is for lower and upper limit of Q to do the intensity summation for Azimuthal plot

- Plot mode: *single*, *multiple* and *waterfall*

Single: 1D plot with slider menu in one hdf5

Multiple: 1D plots **with limited slider menu (only latest opened file)** from multiple hdf5 file

Waterfall: 1D plots in waterfall format within one hdf5 with user defined start, stop, sequence step and plotting gap.


4. Data view window:

Default 1D profile window with multiple choice defined in **Panel 3**

5. Slider menu


Use to cut through the 3D data file in **single plot mode** and through the **latest opened file in multiple file**.

Q-I 1D plot



Browsing through the frames of 1D

Q-I multiple Plots in a single file



Step:


1. Select the a single hdf5 file

2. Select the plotting mode

3. Scroll to the number to be plotted

4. Change legend by right clicking the mouse and choosing "Rename curve" before scrolling to the next curve in panel 3


Q-I multiple Plots from multiple files



Step:

1. Select the plotting mode
2. Select the multiple hdf5 files one by one
3. For each file, scroll to the curve number to be plotted
4. If multiple curves from one hdf5 will be plotted, then legends need to be modified for those curves. Otherwise, no need to change legend between hdf5.

Q-I Waterfall plot from a single file




For 2D “cake” input:
Green shadowed area
is used to define the
integrated azimuthal
angle region

For 1D input:
Green shadowed area
is **NOT** effective

Two modes are available in **log scale** and **linear scale** stacking


Azimuthal plot (only valid for 2D “cake” cut)



Shaded area is used to define q range to be integrated for Azimuthal Intensity

The plotting modes (single, multiple, waterfall (log) and (linear)) are the same to Q-I plot

Simple fitting



1. Click on fit icon to activate the fitting panel.

2. Left click to activate the curve to be fitted


3. Select fitting function

4. Estimate (multiple peaks may be generated automatically)


5. Start fit. The fitting results will be automatically updated

6. Finishing fitting by click dismiss

Plot stacked 1D data into 2D image



Plot 3D data as stacked 2D (“cake” cut) image



Line cut in 2D plot

Choose the line profile style

Integration region (in pixels)

Basic functions are the same as in 1D plot window, e.g. legend, fitting, saving data, etc.

